

OCULUS
Pentacam®
Pentacam® HR
Pentacam® AXL
Anterior Segment
Tomography



OCULUS Pentacam® / Pentacam® HR / Pentacam® AXL

The gold standard for anterior segment tomography



Since its introduction in 2002 the OCULUS Pentacam® has proven to be an indispensable tool for eye care professionals committed to providing precise diagnostics and successful treatment.

Efficient preliminary and follow-up examinations

In a single delegable step, the OCULUS Pentacam® measures the entire anterior eye segment independent of tear film. From the high-resolution Scheimpflug images it calculates a motion-corrected 3D model.

The Pentacam® AXL additionally performs axial length measurements with contact-free optical biometry from the corneal surface to the retina.

Comprehensive analysis

As the basis for precise keratometry, a prerequisite for IOL calculation, and for detection of corneal irregularities as well as wavefront calculation, the Pentacam® provides a complete description of the entire cornea (pachymetry as well as elevation and curvature data).

Scientifically reliable screenings

To perform quick and reliable anterior chamber angle screenings, the OCULUS Pentacam® determines anterior chamber volume, angle and depth automatically and then compares the values in the Fast Screening Report with a normal and a pathological population.

Reliable values for your diagnoses

Using blue light, the OCULUS Pentacam® makes opacities of the cornea, the crystalline lens and IOLs visible so that corneal diseases are detected reliably and cataract progression can be assessed objectively.

The Pentacam® HR stands out by virtue of its bright optics, which make for extremely sharp, highest-quality Scheimpflug images.

Testimonials of satisfied international Pentacam® users

"I cannot imagine any cornea or refractive surgery without a pre and post-op Pentacam® examination."

Dr. Paolo Vinciguerra, MD, Italy



"I have been using the Pentacam® since its inception. I could not imagine taking a refractive patient into surgery without a full tomographic examination. It is the single most important piece of diagnostic equipment in our laser suite. I consider it an indispensable part of my practice."

Dr. Michael W. Belin, MD, USA



"Since the first day I found the Pentacam® to be very reliable and essential for screening refractive candidates for ectasia risk. The Pentacam® is also very important for evaluating patients for corneal diseases, cataract and glaucoma."

Dr. Renato Ambrósio Jr, MD, PhD, Brazil



"In my opinion, the Pentacam® HR is a very useful device for screening before cataract surgery especially for premium IOLs. A detailed assessment can be performed to select the optimal IOL for the respective patient. Especially for toric IOL implantation, the influence of the posterior cornea can be assessed with the Pentacam®."

Dr. Naoyuki Maeda, MD, Japan



"Curvature measurement precision with Pentacam® HR has proved to be better than other tomographers and autokeratometers."

Dr. Jaime Aramberri, MD, Spain



*"The Pentacam® is a must for anyone performing cataract or refractive surgery. No other instrument provides more clinically relevant information." **

Dr. Tobias Neuhann, MD, Germany



* Translated from German by OCLUS.

The Register

Automatic measurement activation, Functions at a glance

Fast Screening Report, General Overview, Topography, Compare Exams

Belin/Ambrósio Display, Corneal staging and progression, Corneal Optical Densitometry, Fourier Analysis, Corneal Rings

Cataract Pre-OP Display, Corneal Power Distribution, Zernike Analysis, PNS and 3D Cataract Analysis, IOL power calculation

CSP Report, 3D pIOL Simulation and Aging Prediction, Contact Lens Fitting

Software overview, Networkability

Basic principles

Screening

Refractive

Cataract

More software

Configuration & Tech.


Anterior Segment Tomography

Fast, reproducible, delegable

Thanks to automatic measurement activation the Pentacam® gives you an overall view of the anterior eye segment in two seconds. Measurements are made independent of tear film and examiner. The Pentacam® HR evaluates up to 138 000 measuring points.

The most important functions at a glance:

- Topography of anterior and posterior corneal surface
- Full-surface pachymetry
- Evaluating the corneal shape
- 3D anterior chamber analysis
- Fast Screening Report with data for crucial parameters
- Normative corneal wavefront data
- Total Corneal Refractive Power (TCRP)
- 3D densitometry of the cornea and the crystalline lens
- Contact lens fitting software
- Anterior segment tomography
- Measurement of the cornea scleral profile (CSP)
- Zernike and Fourier Analysis
- Axial length measurements (Pentacam® AXL only)



The Pentacam®: Indispensable for you and your patients

Before the surgery, my doctor explained what the Pentacam® examination involved. I was amazed that I was still even able to see at all with such a clouded lens!



Cataract

Benefit from simple, comprehensive corneal screenings. Fast Screening Report, early detection of corneal abnormalities according to Belin/Ambrósio and evaluation of corneal optical densitometry. The Cataract Pre-OP Display provides an overview of all criteria for reliable selection of a premium IOL. The IOL Calculator allows calculating IOLs for treated as well as untreated corneas, giving due consideration to posterior surface measurements and total corneal refractive power. Prof. Olsen's Ray-Tracing formula also allows calculating difficult cases, such as with irregular corneal geometry.

I simply didn't want to wear spectacles anymore! During the pre-op examination the Pentacam® results gave me certainty.



Refractive assessment

The Belin/Ambrósio Display supports you in early detection of corneal abnormalities. A final parameter is calculated and represented in color. This is how the Pentacam® software assists you in making diagnoses. Keratometry, asphericity and full-surface pachymetry assist you in planning refractive surgery. The Pentacam® measures irregular corneas with great precision, presenting the necessary parameters in a clear and application-oriented fashion. This ensures reliable planning for implantation of corneal rings, CXL and cornea transplantations.

My risk of developing glaucoma was detected by my ophthalmologist early on. He was able to initiate appropriate measures immediately.



Glaucoma

Glaucoma is one of the most common eye diseases. Make use of the evaluation of anterior chamber angle and volume based on published papers and clinical pictures. In particular, studies* have identified automatically calculated anterior chamber volume as a sensitive parameter.

* Assessment of the anterior chamber parameters after laser iridotomy in primary angle close suspect using Pentacam® and gonioscopy; Alireza et al; Int J Ophthalmol, 2013, 6(5):680-684
Comparison of scheinpflug imaging and spectral domain anterior segment optical coherence tomography for detection of narrow anterior chamber angles; Grewal et al; Eye Vol.: 25

Basic principles

Screening

Refractive

Cataract

More software

Configuration & Tech.

General Screening

Filter, evaluate, represent

The great challenge in clinical routine is filtering, evaluating, and representing data in a clear-cut way. This is exactly what the Fast Screening Report does. The crucial data are represented in such a way that you can gain a comprehensive picture of your patient at a glance.



Fast Screening Report

- Data are gathered from published studies and stored in the Pentacam® software.
- The distribution of normal values in a population is represented by the bars shaded in gray.
- The diagrams show the distribution for normal (green) and pathological (red) eyes.
- Sources which the evaluations are based on are cited to provide additional information.
- In cases of irregularities the individual displays for detailed findings appear in the interactive navigation bar.

Early detection of corneal abnormalities

Individual patient value


Distribution in normal population

Distribution of a normal (green) and a pathological (red) population

Interactive navigation bar


Screening and follow-up exams ensure safety

After initial information has been provided by the Fast Screening Report additional displays show further, case-specific information. Depending on the irregularities the interactive navigation bar recommends displays which lead you to suitable detailed analyses – custom-tailored to the individual patient.




General Overview

The Scheimpflug images allow you to make a qualitative assessment of the anterior eye segment. Opacities of the cornea or the crystalline lens are represented impressively. Keratometry, pachymetry and asphericity allow for an initial analysis of the corneal surface. The anterior chamber is defined by chamber volume, angle and depth. Intraocular pressure is corrected on the basis of central corneal thickness.



Compare 2 Exams

The comparative displays in the Pentacam® software allow following up on screenings, performing pre-op to post-op differential analyses and illustrating one's explanations to patients. All Pentacam® color maps are available for comparison. There is one display for comparing two and another for comparing four examinations.



4 Maps Refractive

Topography and elevation maps make quantitative assessments of the corneal surface possible. Especially the posterior elevation map shows pathological changes very early on. The corneal thickness map clearly shows the thinnest position to help you plan refractive surgery, for example. Combining elevation data of the anterior and posterior corneal surface, topography and pachymetry will help you detect abnormalities early on.

Screening the Corneal Shape


Early detection of corneal abnormalities

The Pentacam® is equipped with intuitive and user-friendly software features to ensure patient safety and help physicians select the best options for optimal results.

Belin/Ambrósio Display

This screening is geared towards early detection of corneal irregularities. Using this program, structural pachymetric progression and the anterior and posterior corneal surface are evaluated. It also has a database for myopic and hyperopic eyes. All the individual parameters are consolidated in one final factor based on regression analysis.

The individual parameters and the final factor are represented in the colors white, yellow, and red to indicate the degree of abnormality.



Calculation and analysis of standard and enhanced anterior elevation maps

Calculation and analysis of standard and enhanced posterior elevation maps



Deviation of individual parameters from normal data

Pachymetric progression

Final factor for early detection of corneal irregularities

Corneal staging and progression assessment

The Pentacam® assists you with state-of-the-art software modules for efficient corneal staging and progression.



Topometric/KC Staging

The Topometric/KC Staging display is used for assessing corneal irregularities. It contains the Belin ABCD Corneal Staging system, a classification method developed by Prof. Michael Belin, USA¹. It shows all parameters relevant to corneal-shape classification at a glance, including the posterior corneal surface and the thinnest pachymetric values. This staging system can also be used retrospectively for classification of a patient's earlier examination data.

1 Global consensus on keratoconus and ectatic diseases; Gomes JA1, Tan D, Rapuano CJ, Belin MW, Ambrósio R Jr, Guell JL, Malecaze F, Nishida K, Sangwan VS; Group of Panelists for the Global Delphi Panel of Keratoconus and Ectatic Diseases; Cornea. 2015 Apr; 34(4): 359-69.

Belin ABCD Progression Display

The evaluation is based on the Belin ABCD Corneal Staging module and can also be used retrospectively. Four different parameters are represented graphically for easy determination of any progression or regression trend (following crosslinking). Changes found are shown with confidence intervals representing normal or pathological populations.

1 Large Color Map

This display can be used, amongst other purposes, to separately view the posterior elevation map calculated from the internal 3D model of the Pentacam® software. This is made possible by superimposing a reference surface onto the posterior corneal surface, bringing its irregularities into clearer view.


Refractive Screening

Planning refractive laser surgery

Corneal assessment involves more than mere topography. It rather calls for a holistic evaluation. The Pentacam® represents pachymetric progression, allowing for an evaluation of corneal structure. Corneal optical densitometry facilitates targeted slit lamp examinations that enable you to detect diseases early on.

Corneal Optical Densitometry

Panorama images of the cornea visualize corneal irregularities. This allows for objective quantification and follow-up. Optical densitometry can be evaluated using a table or a color chart. In the table the measured values are represented according to various zones and layers. Corneal optical densitometry results are displayed in reference to published age-related normative data. This provides a basis for the detection of other diseases such as Fuchs' Dystrophy.



Corneal densitometry of selected layer

Horizontal Scheimpflug image

Vertical Scheimpflug image

Graphs for average and maximum densitometry of selected layer




Average or maximum densitometry map

Graphs for average and maximum densitometry of selected ring segment

Densitometry values of different layers and ring segments

Assess the entire cornea qualitatively and quantitatively

The Pentacam® measures the entire cornea from limbus to limbus non-invasively. This ingenious method is unaffected by the tear film and guarantees accurate measurements of the corneal center. It provides a structured representation of all the parameters crucial for planning refractive surgery.



Refractive

This display provides a neat overview of all parameters important to the refractive surgeon. It shows pachymetry progression as well as anterior corneal curvature and eccentricity values.

Fourier Analysis

The refractive power of the front surface of the cornea consists of different components. The Fourier Analysis identifies four of them:

- Spherical component
- Decentration
- Regular astigmatism
- Irregularities

Pathological changes can be quantified and possible effects on visual acuity can be explained.

Corneal Rings

This display shows all the parameters necessary for planning corneal ring implantation. Dependent on the selected surgical procedure – manual dissection technique or femtosecond laser – corneal thickness is represented in specific areas and segments.

Easily Find the Perfect IOL, Premium IOL, or ICL

Fine-tune your patient's vision


The Pentacam® measures TCRP with ray tracing, giving due consideration to the front and back corneal surfaces. The Cataract Pre-OP Display summarizes all data relevant to your decisions in selecting toric, aspheric, multifocal, or standard IOLs for your patients.

Cataract Pre-OP Display

The Cataract Pre-OP Display was developed in collaboration with Prof. Dr. Naoyuki Maeda from the University Medical School in Osaka, Japan.

This display assists users in:

- Cataract pre-op screening
- IOL selection incl. premium IOL
- IOL power calculation
- Determining criteria for ICL selection



Preoperative diagnostics is your success concept

For documentation of cataract progression and optimal surgery planning the Pentacam® provides cataract surgeons with comprehensive analysis options.




Corneal Power Distribution

The table shows the refractive power of the cornea in various zones and rings. This allows for individual assessment of the posterior corneal surface influence on the total corneal power, astigmatism and axis.



Zernike Analysis

The wavefront error of the cornea as a whole is calculated individually using ray tracing. Higher order aberrations are calculated and represented in comparison to a normal population.



PNS and 3D Cataract Analysis

Through blue-light illumination opacities of the natural lens become visible. Two- and three-dimensional quantification of lens opacification in the individual layers and of posterior capsular opacification is performed. The PNS (Pentacam® Nucleus Staging) function permits optimized ultrasonic energy output in phacoemulsification as well as optimized effective phaco time (EPT) in Femto assisted cataract surgery.

IOL Power Calculation

IOL power calculation for virgin and post corneal refractive eyes

The IOL Calculator provides different formulas for virgin and post corneal refractive eyes. Spheric, aspheric, multifocal and toric IOLs are included in the database. The data of the implanted IOL as well as the post-op refraction can be stored with only two clicks. The software also includes IOL constant optimization (Prof. Wolfgang Haigis).

TCRP for toric IOL power calculation

Toric IOLs can be calculated based on the individual SIA (surgically induced astigmatism) and TCRP. The latter takes the posterior corneal surface into account, improving the accuracy of IOL power calculation¹. Included formulas: Hoffer Q, Holladay I, SRK/T, Haigis, Barrett Universal II, Savini¹ and Barrett for toric IOLs, HillPotvinShammas² after myopic lasik, HillPotvin after RK³, double K of Hoffer, SRK/T and Holladay and the Olsen formula (ray tracing). The Holladay consultant and OKULIX (ray tracing) can be linked to the software.

- 1 An Analysis of the Factors Influencing the Residual Refractive Astigmatism After Cataract Surgery With Toric Intraocular Lenses; Giacomo Savini and Kristian Naeser; IOVS2015; 56:827-835.DOI:10.1167/iov.14-15903
- 2 New algorithm for intraocular lens power calculations after myopic laser in situ keratomileusis based on rotating Scheimpflug camera data; Richard Potvin, OD, Warren Hill, MD; JCRS 2015; 41:339-347 Q 2015 ASCRS and ESCRS
- 3 New algorithm for post-radial keratotomy intraocular lens power calculations based on rotating Scheimpflug camera data; Richard Potvin, OD, Warren Hill, MD; JCRS 2013; 39:358-365 Q 2012 ASCRS and ESCRS

! This software module is available for the Pentacam® HR (optional) and the Pentacam® AXL only



Axial length

SIA and target refraction SEQ



Automatic plausibility check

Implanted and stored IOL and surgery data



IOL power calculation
Post-op refraction data for IOL constant optimization

IOL power calculation printouts

All necessary information, such as IOL power results and resulting SEQ, is summarized in a clear and intuitive printout, ready for use right before surgery.

OCULUS 			
Patient 01 Demo Cataract, Cataracta Centralis	Exam Date 17.10.2006		
D.o.B. 06.04.1935	Exam Time 11:20:06		
ID OD grade 2, OS opaque epithel	Planning Date 25.03.2019		
Surgeon TestSurgeon, Test	Eye Status		
Cor./Diseases: Info:			
		OS (Left)	
QS OK	Tgt Refr SEQ -0.25 D	SIA 0.15 D @ 90°	
AL (opt.) manual 23.600 mm (M)	SimK 15° (n = 1.3375)		
ACD (Ext.) 3.73 mm	K1 43.1 D / 7.83 mm @ 149°		
Pupil Dia 2.65 mm	K2 43.4 D / 7.78 mm @ 59°		
HWTW	K Avg 43.2 D / 7.80 mm		
Chord μ 0.11 mm	Astig 0.3 D		
TCRP WFAZ40, 6mm 0.393 μm	K1 Pre-Refr.-Surg.		
TCRP WFAHOA, 4mm 0.097 μm	K2 Pre-Refr.-Surg.		
OS-OD Test OK			
1 Alcon AcrySof IQ SN60WF		2 Alcon AcrySof IQ SN60WF	
Haigis K1/K2 (SimK 15°): 43.2 D (nRef=1.3375) IOL SEQ Emm. = +21.76 D a0:-0.769 a1:0.234 a2:0.217		Barrett Universal II K1/K2 (SimK 15°): 43.2 D IOL SEQ Emm. = +21.59 D ABarrett: 119	
IOL SEQ	Refraction SEQ	IOL SEQ	Refraction SEQ
+21.50	+0.18	+21.00	+0.41
+22.00	-0.17	+21.50	+0.06
+22.50	-0.52	+22.00	-0.29
+23.00	-0.87	+22.50	-0.65
+23.50	-1.23	+23.00	-1.01
3 Alcon AcrySof IQ SN60WF		4 Alcon AcrySof IQ SN60WF	
Holladay1 K1/K2 (SimK 15°): 43.2 D IOL SEQ Emm. = +21.34 D SF: 1.84		SRK/T K1/K2 (SimK 15°): 43.2 D IOL SEQ Emm. = +21.25 D ASRKT: 119	
IOL SEQ	Refraction SEQ	IOL SEQ	Refraction SEQ
+21.00	+0.23	+21.00	+0.17
+21.50	-0.11	+21.50	-0.17
+22.00	-0.45	+22.00	-0.52
+22.50	-0.79	+22.50	-0.87
+23.00	-1.14	+23.00	-1.22

Printout of an IOL power calculation for an aspherical IOL

OCULUS 			
Patient 01 Demo Cataract, Cataracta Centralis	Exam Date 17.10.2006		
D.o.B. 06.04.1935	Exam Time 11:20:06		
ID OD grade 2, OS opaque epithel	Planning Date 25.03.2019		
Surgeon TestSurgeon, Test	Eye Status		
Cor./Diseases: Info:			
		OS (Left)	
QS OK	Tgt Refr SEQ -0.25 D	SIA 0.15 D @ 90°	
AL (opt.) manual 23.600 mm (M)	SimK 15° (n = 1.3375)		
ACD (Ext.) 3.73 mm	K1 43.1 D / 7.83 mm @ 149°	TCRP 3mm, zone, pup. 42.6 D @ 135°	
Pupil Dia 2.65 mm	K2 43.4 D / 7.78 mm @ 59°	42.7 D @ 45°	
HWTW	K Avg 43.2 D / 7.80 mm	42.7 D	
Chord μ 0.11 mm	Astig 0.3 D	0.1 D	
TCRP WFAZ40, 6mm 0.393 μm	K1 Pre-Refr.-Surg.		
TCRP WFAHOA, 4mm 0.097 μm	K2 Pre-Refr.-Surg.		
OS-OD Test OK			
13 Alcon AcrySof IQ Toric SN6AT(2-9)			
Haigis K1/K2 (SimK 15°): KAvg = 43.2 D (n=1.3375) IOL SEQ Emm. = +22.08 D a0:-0.323 a1:0.213 a2:0.208		Savini Toric TCRP 3mm, zone, pup.: Astig: 0.1D @ 45° IOL Toricity Emm. = +0.26 D @ 17°	
IOL SEQ	Refraction SEQ	IOL Toricity	Astig. Res.
+21.50	+0.40	n.a.	
+22.00	+0.05	n.a.	
+22.50	-0.29	T2 1.00 D	+0.50 D @ 107°
+23.00	-0.64	T3 1.50 D	+0.84 D @ 107°
+23.50	-0.99	T4 2.25 D	+1.35 D @ 107°
IOL	SEQ: +22.50 D T2 (1.00 D)		
IOL Axis	17°		
Residual Refr.	SEQ -0.29 C +0.50 @ 107°		
Incision Axis	90°		

Printout of an IOL power calculation for a toric IOL

IOL printout

For optimal pre-op decision making, all relevant data (i.a. ACD, HWTW) is displayed, the following of which deserve particular mention:

- Chord μ
- Total WFA Z40, 6 mm
- Total WFA HOA, 4 mm

The result of the automatic plausibility check is displayed, corroborating or querying the measured data to support you in your decision making.

! This software module is available for the Pentacam® HR (optional) and the Pentacam® AXL only

Basic principles

Screening

Refractive


Cataract

More software

Configuration & Tech.


More Software Modules

Various software options for different challenges



Contact Lens Fitting


The Pentacam® provides tear film independent measurements, offering optimal support in fitting contact lenses. Its dynamic fluo image simulation makes it possible to predict the seat as well as any inclination or displacement of the contact lens. Up to seven simulations can be displayed simultaneously.



CSP Report

The cornea scleral profile (CSP) is measured with a coverage of up to 18 mm, capturing not only the cornea but also the sclera. Together with the sagittal height data this provides you with parameters that are important for fitting scleral lenses. All 250 Scheimpflug images are generated tear film independently and without the patients having to move their eyes. Thanks to user-independent image release the CSP Report is just as reproducible as any other measurement taken with the Pentacam®.

! This software module is available for the Pentacam® and the Pentacam® HR only



3D pIOL Simulation and Aging Prediction

This software module assists you in preoperative planning of iris-fixated phakic IOLs. After the subjective refraction is entered the software calculates the necessary pIOL refractive power dependent on the selected pIOL type. The position of the pIOL in the anterior chamber is automatically calculated in 3D and represented in the Scheimpflug images. The minimal distances between the pIOL and the crystalline lens as well as the endothelium are calculated automatically in 3D and displayed numerically as well as in a color map.

! This software module is available for the Pentacam® HR (optional) and the Pentacam® AXL only

All Features at a Glance

Customize the OCULUS Pentacam® / Pentacam® HR / Pentacam® AXL to your own requirement

Software included	Pentacam®	Pentacam® HR	Pentacam® AXL
General Overview	✓	✓	✓
Fast Screening Report	✓	✓	✓
Topometric / KC-Staging (Belin ABCD Corneal Staging)	✓	✓	✓
Belin/Ambrósio Display	✓	✓	✓
Belin ABCD Progression Display	✓	✓	✓
Elevations and Topography of corneal front and back surface	✓	✓	✓
4 Maps Refractive	✓	✓	✓
Iris Image and HWTW	✓	✓	✓
Compare 2 Exams	✓	✓	✓
Tomography	✓	✓	✓
Scheimpflug Image Overview	✓	✓	✓
3D anterior chamber analysis	✓	✓	✓
Compare 2 Exams Scheimpflug Images	✓	✓	✓
Contact Lens Fitting incl. Fourier Analysis	✓	✓	✓
Additional software	Pentacam®	Pentacam® HR	Pentacam® AXL
Holladay Report	○	✓	✓
CSP Report	○	○	-
3D pIOL Simulation and Aging Prediction	-	○	○
IOL Calculator	-	○	✓
Corneal Optical Densitometry	○	✓	✓
PNS and 3D Cataract Analysis	○	✓	✓

✓ included ○ optional - not available

Optional examination functions (included with the Pentacam® AXL)

Cataract Software Package

- Cataract Pre-OP Display
- Corneal Power Distribution
- Zernike Analysis with normative corneal wavefront
- Keratometry overlay in Iris Image
- PNS and 3D Cataract Analysis
- Total Corneal Refractive Power (TCRP)
- Automatic calculation of the anterior chamber angle in 360°, measurement based on Scheimpflug images
- 4 Maps Topometric and 4 Maps Chamber
- Show 2 Exams
- Compare 4 Exams
- True Net Power (TNP)

Refractive Software Package

- Corneal Optical Densitometry
- Corneal Rings
- Fourier Analysis
- Refractive
- Pachymetric
- 4 Maps Selectable
- Show 2 Exams
- Compare 4 Exams

Optional evaluation functions

DICOM Interface

Basic principles

Screening

Refractive

Cataract

More software

Configuration & Tech.

Scheimpflug Images

Impressive, convincing, simply sharp



Descemet Stripping After Endothelium Keratoplasty (DSAEK)




Extremely thin keratoconus cornea




Penetrating Keratoplasty (PKP)



Implantable Collamer® Lens (ICL)



Nuclear sclerosis



Narrow (anterior) chamber angle

The Fascination of Technology


Ingenious yet simple

The Floating License Key – for maximum flexibility

This basic Pentacam® software is already available at all workplaces in your network. You decide which optional examination and evaluation functions you need in addition. You also choose the number of optional evaluation software features which should be available simultaneously. The Floating License Key (FLK) activates the corresponding licenses and makes them available on your network. To help you decide what you need all the optional evaluation functions can be accessed 20 times for demonstration purposes.

Efficiency and productivity through networking

The OCULUS Patient Data Management system (PDM) optimizes your work processes. It is always included in the scope of delivery; it organizes patient and examination data from all OCULUS instruments. The PDM is network-compatible and can be incorporated into many Electronic Medical Record (EMR) systems. Needless to say, the OCULUS PDM communicates with the DICOM environment and makes results available in DICOM format.



Basic principles

Screening

Refractive

Cataract

More software

Configuration & Tech.

Pentacam®/Pentacam® HR /Pentacam® AXL

Technical Data

Scheimpflug camera	Pentacam®	Pentacam® HR	Pentacam® AXL
Camera	digital CCD camera	digital CCD camera	digital CCD camera
Light source	blue LED (475 nm UV-free)	blue LED (475 nm UV-free)	blue LED (475 nm UV-free)
Processor	DSP with 400m operations/s	DSP with 400m operations/s	DSP with 400m operations/s
Speed	50 images in 2 seconds ¹⁾	100 images in 2 seconds ²⁾	100 images in 2 seconds ²⁾
Measurement range	Pentacam®	Pentacam® HR	Pentacam® AXL
Curvature	3 - 38 mm 9 - 99 D	3 - 38 mm 9 - 99 D	3 - 38 mm 9 - 99 D
Precision	± 0.2 D	± 0.1 D	± 0.1 D
Reproducibility	± 0.2 D	± 0.1 D	± 0.1 D
Operating distance	80 mm (3.1 in)	80 mm (3.1 in)	80 mm (3.1 in)
Axial length	-	-	14 - 40 mm
Technical specifications	Pentacam®	Pentacam® HR	Pentacam® AXL
Dimensions (W x D x H)	275 x 320 - 400 x 500 - 530 mm (10.8 x 12.6 - 15.7 x 19.7 - 20.9 in)	275 x 320 - 400 x 500 - 530 mm (10.8 x 12.6 - 15.7 x 19.7 - 20.9 in)	278 x 320 - 400 x 500 - 532 mm (10.8 x 12.6 - 15.7 x 19.7 - 20.9 in)
Weight (measuring head)	7.2 kg	7.8 kg	8.4 kg
Max. power consumption	35 W	42 W	42 W
Recommended computer specifications	Intel® Core™ i5, 500 GB HDD, 8 GB RAM, Windows® 10, Intel® HD Graphics DIN EN 62368-1 or DIN EN 60950		

¹⁾ Scheimpflug image of the entire anterior segment

²⁾ Cornea Fine Scan

CE CE in accordance with Regulation (EU) 2017/745 on Medical Devices



WWW.OCULUSUSA.COM




OCULUS is certified by TÜV according to
DIN EN ISO 13485 MDSAP

Distributed in Canada By:
Clarion Medical Technologies
125 Fleming Drive, Cambridge, Ontario, N1T 2B8
Tel. +1 519-620-3900 • +1 800-668-5236
Email: info@clarionmedical.com
Website: www.clarionmedical.com

Headquarters:
OCULUS Optikgeräte GmbH
Postfach • 35549 Wetzlar • GERMANY
Tel. +49 641 2005-0 • Fax +49 641 2005-295
Email: export@oculus.de • www.oculus.de

OCULUS USA:
OCULUS, Inc.
17721 59th Avenue NE • Arlington • WA 98223
Tel. +1 425 670 9977 • Fax +1 425 670 0742
Email: sales@oculususa.com • www.oculususa.com



US version. The availability of products and features may vary by country. OCULUS reserves the right to change product specifications and design. All information is valid at the time of printing (05/21).