



ENDYMED™ | PURE


More Powerful, More Possibilities

2.0


The PURE 2.0, is the latest advance in the evolution of RF treatment systems, and is powered by unique, FDA cleared 3DEEP RF technology. Unlike Monopolar or Bipolar RF technology, which is more superficial and can often be painful or risky to the epidermis, 3DEEP's delivers controlled energy to the depth of the dermis, reaching 52-55C, with minimal epidermal heat. This is achieved through gently focusing the device energy into the dermis for optimal, painless collagen remodelling.

EnyMed's 3DEEP technology incorporates real time skin impedance monitoring, providing the most optimal energy to the treatment target area by making continuous adjustments based on the tissue impedance over the current treatment area.


3DEEP - Phased Controlled RF Technology



A **Bipolar RF**
Relatively superficial
and less efficient



B **Monopolar RF**
Energy flows
uncontrolled in the
tissues, requiring
intense cooling



C **3DEEP RF**
Focused and contained
deep energy flow with
no need for cooling

ENDYMED™ | iFine iFINE HANDPIECE:

The small head enables you to treat comfortably within the periorbital area, while the ergonomic design results in less user fatigue.

